Keron

FINE QUALITY PRODUCTS

Nitrie Profi

REF	size
15380	6,5- 7/5
15381	7,5- 8/ M
15382	8,5- 9/ L
15383	9,5- 10/ XL
15384	10,5- 11/ XXL

DE Kategorie III - Chemische und biologische Risiken -
FR Catégorie III - Risques chimiques et biologiques -
EN Category III - Chemical and biological risks -
IT Categoria III - Rischi chimici e biologici -
PT Categoria III - Riscos químicos e biológicos -
ES Categoría III - Riesgos químicos e biológicos -
NL Kategorie III - Chemische en biologische risico's -
SV Kategori III - Kemiska och biologiska risker -
FI Luokkia III - Kemialliset ja biologist vaarat -
DA Kategori III - Kemiske og biologiske risici -
NO Kategori III - Kjemiske og biologiske risikoer -
PL Kategorie III - Zagrozenia chemiczne i biologiczne -
HU III. kategória - Kémiai és biológiai kockázatok -
SK Kategória III - Chemické a biologické riziká -
CS Kategorie III - Chemická a biologická rizika -
SL Kategorija III - Kemična in biološka tveganja -
HR Kategorija III - Kemijske i biološke opasnosti -
RU Категория III - Химические и биологические риски -
RO Categoria III - Kimyasal ve biyolojik riskler -
TR Kategorie III - Riskli biyolojik ve kimyasal

DE Kategorie III - Chemische und biologische Risiken -
FR Catégorie III - Risques chimiques et biologiques -
EN Category III - Chemical and biological risks -
IT Categoria III - Rischi chimici e biologici -
PT Categoria III - Riscos químicos e biológicos -
ES Categoría III - Riesgos químicos e biológicos -
NL Kategorie III - Chemische en biologische risico's -
SV Kategori III - Kemiska och biologiska risker -
FI Luokkia III - Kemialliset ja biologist vaarat -
DA Kategori III - Kemiske og biologiske risici -
NO Kategori III - Kjemiske og biologiske risikoer -
PL Kategorie III - Zagrozenia chemiczne i biologiczne -
HU III. kategória - Kémiai és biológiai kockázatok -
SK Kategória III - Chemické a biologické riziká -
CS Kategorie III - Chemická a biologická rizika -
SL Kategorija III - Kemična in biološka tveganja -
HR Kategorija III - Kemijske i biološke opasnosti -
RU Категория III - Химические и биологические риски -
RO Categoria III - Kimyasal ve biyolojik riskler -
TR Kategorie III - Riskli biyolojik ve kimyasal

EN ISO 374-5:2016	EN ISO 374-1:2016 +A1:2018 Type A
	
VIRUS	JKOPST

DE PPE unterliegt dem Konformitätsbewertungsverfahren Modul C2 unter Überwachung der benannten Stelle -
FR L'EPI est soumis à la procédure d'évaluation de conformité Module C2, sous surveillance de l'organisme notifié -
EN PPE is subject to the conformity assessment procedure Module C2 under surveillance of the notified body -
IT Il DPI è soggetto alla procedura di valutazione della conformità Modulo C2 sotto la sorveglianza dell'organismo accreditato -
NL PPE is onderhevig aan de procedure voor conformiteitsbeoordeling Module C2, onder toezicht van de aangemelde instantie -
ES El EPI está sujeto al procedimiento de evaluación de la conformidad (módulo C2) bajo la supervisión del organismo notificado -
PT PPE sujeito ao procedimento de avaliação de conformidade de, módulo C2, sob a supervisão do organismo notificado -
SV Personell skyddsutrustning omfattas av förfarandet för bedömning av överensstämmelse, modul C2 under övervakning av det anmälda organet -
FI Henkilösuojaimiin (PPE) kohdistetaan yhdennämuksisuusarviointin moduulin C2 mukaisesti ilmoitetun laitoksen valvonnassa -
DA Det personlige værnemiddel er omfattet af proceduren for overensstemmelsesvurdering modul C2 under overvågning af det bemyndede organ -
NO PPU er ikke underlagt prosedyre for vurdering av samsvar modul C2 ved overvåkning av teknisk kontrollorgan -
PL Wyposażenie ochrony osobistej podlega procedurze oceny zgodności wg modułu C2 pod nadzorem jednostki notyfikowanej -
HR EPI PPE (personalni zaštitni oprema - egipveni védőeszköz) a C2 modul megfelelőségértékelési eljárásának hatálya alá tartozik a bejelentett szervezet irányításával -
SK PPE podlieha postupu posudzovania zhody Moduľ C2 pod dohľadom notifikovaného orgánu -
CS Osobní ochranné prostředky podléhají postupu posuzování shody modulu C2 pod dohledem notifikovaného orgánu -
SL Za OVO velja modul C2 v postopku ugotavljanja skladnosti pod nadzorom prijavljenega organa -
HR PPE je predmet postupka procjene skladnosti Module C2 pod nadzorom certifikacijskog tijela -
RU СИЗ проходит процедуру оценки соответствия по модулю C2 под надзором официального органа -
RO EPP este supus Modulului C2 al procedurii de evaluare a conformității sub supravegherea organismului notificat -
TR PPE, onaylanmış kurumun gözetimli altında uygunluk değerlendirmeye prosedürü Modül C2'ye tabidir.

SATRA Technology Europe Ltd
Braceown Business Park
Crazeon, D15YN2P, Ireland
Notified Body No.: 2777

DE Baumusterprüfung durch: -
FR Examen de type effectué par: -
EN EU type examination carried out by: -
IT Esame del tipo mediante: -
ES Examen de tipo por: -
PT Exame UE de tipo realizado por: -
NL Typeonderzoek door: -
SV Typkontroll genom: -
FI Tyypitarkastus: -
DA Typeafprøvning af: -
NO Typeprøving gjennom: -
PL Badanie typu przez: -
HU A típusvizsgálatot végzte: -
SK Skúška konštrukčného zvoru prostredníctvom: -
CS Přezkoušení typu prostřednictvím: -
SL Testiranje vzorca opravljeno pri: -
HR Tipsko ispitivanje po: -
RU Испытание опытного образца: -
RO Examinare de tip de către: -
TR Yapi numunesi testineşi gerçekleştirilen:

EU: SATRA Technology Europe Ltd
Braceown Business Park
Crazeon, D15YN2P, Ireland
Notified Body No.: 2777

SATRA Technology Europe Ltd
Braceown Business Park
Crazeon, D15YN2P, Ireland
Notified Body No.: 2777

DE Baumusterprüfung durch: -
FR Examen de type effectué par: -
EN EU type examination carried out by: -
IT Esame del tipo mediante: -
ES Examen de tipo por: -
PT Exame UE de tipo realizado por: -
NL Typeonderzoek door: -
SV Typkontroll genom: -
FI Tyypitarkastus: -
DA Typeafprøvning af: -
NO Typeprøving gjennom: -
PL Badanie typu przez: -
HU A típusvizsgálatot végzte: -
SK Skúška konštrukčného zvoru prostredníctvom: -
CS Přezkoušení typu prostřednictvím: -
SL Testiranje vzorca opravljeno pri: -
HR Tipsko ispitivanje po: -
RU Испытание опытного образца: -
RO Examinare de tip de către: -
TR Yapi numunesi testineşi gerçekleştirilen:

DE Baumusterprüfung durch: -
FR Examen de type effectué par: -
EN EU type examination carried out by: -
IT Esame del tipo mediante: -
ES Examen de tipo por: -
PT Exame UE de tipo realizado por: -
NL Typeonderzoek door: -
SV Typkontroll genom: -
FI Tyypitarkastus: -
DA Typeafprøvning af: -
NO Typeprøving gjennom: -
PL Badanie typu przez: -
HU A típusvizsgálatot végzte: -
SK Skúška konštrukčného zvoru prostredníctvom: -
CS Přezkoušení typu prostřednictvím: -
SL Testiranje vzorca opravljeno pri: -
HR Tipsko ispitivanje po: -
RU Испытание опытного образца: -
RO Examinare de tip de către: -
TR Yapi numunesi testineşi gerçekleştirilen:

EU: SATRA Technology Europe Ltd
Braceown Business Park
Crazeon, D15YN2P, Ireland
Notified Body No.: 2777

SATRA Technology Europe Ltd
Braceown Business Park
Crazeon, D15YN2P, Ireland
Notified Body No.: 2777

DE Baumusterprüfung durch: -
FR Examen de type effectué par: -
EN EU type examination carried out by: -
IT Esame del tipo mediante: -
ES Examen de tipo por: -
PT Exame UE de tipo realizado por: -
NL Typeonderzoek door: -
SV Typkontroll genom: -
FI Tyypitarkastus: -
DA Typeafprøvning af: -
NO Typeprøving gjennom: -
PL Badanie typu przez: -
HU A típusvizsgálatot végzte: -
SK Skúška konštrukčného zvoru prostredníctvom: -
CS Přezkoušení typu prostřednictvím: -
SL Testiranje vzorca opravljeno pri: -
HR Tipsko ispitivanje po: -
RU Испытание опытного образца: -
RO Examinare de tip de către: -
TR Yapi numunesi testineşi gerçekleştirilen:

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Ⓢ Instructions for use

The glove types named here meet with the requirements of Personal Protection Equipment Regulation (EU) 2016/425 and Regulation 2016/425 as brought into UK law and amended, regulation 2017/745, UK MDR 2021, and is compliant with the above mentioned harmo-nised designated standards EN4551/-4, EN ISO 374-1:2016+A1:2018 and EN ISO 21420: 2020+A1:2024. As the gloves are designed to special purposes, the lengths may deviate from the requirements EN ISO 21420: 2020+A1:2024.

Results chemical tests :
n-heptane (I): Level 6 30% Hydrogen Peroxide (P): Level 3
40% Sodium Hydroxide (K): Level 6 40% Hydrofluoric Acid (S): Level 2
65% Nitric Acid (M): Level 1 37% Formaldehyde (T): Level 6
25% Ammonium Hydroxide (O) Level 2

Level	1	2	3	4	5	6
Breakthrough time (min)	> 10	> 30	> 60	> 120	> 240	> 480

EN ISO 374-4:2019 - Resistant to degradation against chemicals:
n-heptane (I): 27,30 % 30% Hydrogen Peroxide (P): -11,9 %
40% Sodium Hydroxide (K): -11,5 % 40% Hydrofluoric Acid (S): X
65% Nitric Acid (M): 97,8 % 37% Formaldehyde (T): -3,9 %
25% Ammonium Hydroxide (O) 18,80 %

EN ISO 374-5:2016: Protection against bacteria and fungi: Pass
Protection against viruses: Pass

Warning: The degradation performance for 40% Hydrofluoric Acid is unknown, it is therefore important to carefully monitor the glove's physical characteristics during use and stop use immediately if there are signs of degradation.

Warning / risk assessment: a) Gloves to protect substances and mixtures which are hazardous to health, and harmful biological agents. Important: The gloves are recommended for use in situations where only low chemical protection is identified as needed. b) While selecting an equipment, user should perform risk analysis based on the intended use and determine the suitability based on product's test standards and protection levels obtained. c) Information provide does not reflect the actual duration of protection in the workplace due to other factors influencing the performance, such as temperature, abrasion, and degradation, and the difference between mixtures and pure chemicals d) Information regarding protection refers to the working surface, i.e. the 'palm' of the glove, which has been submitted to testing. e) Gloves should be thoroughly inspected for damages (specially for knits and holes) before use. If any damage is found avoid usage. f) The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. g) It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation h) When used, protective gloves may provide less resistance to the chemical tested due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc, may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves i) The maximum wear time depends on the activity being carried out and the person. j) EN ISO 374-4:2019 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. k) The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. l) This product contains Nitrile Butadiene Rubber and compounding chemicals which may cause allergic reaction to some individuals. m) The glove offers no protection against perforation with sharp objects, e.g. injection needles. **Donning:** Insert the hand into the glove carefully without damaging the glove. **Doffing:** - Grasp the outside of the glove from the wrist area • Peel the glove away from the hand, hold it in the opposite gloved hand • Slide an un-gloved finger under the wrist of the remaining glove, being careful not to touch the contaminated surface of the glove • Additional information will be supplied upon request. The gloves can be used in chemical industry, petrochemical industry, auto motive industry, aircraft industry and facility maintenance. **Storage:** To be stored in cool dry place away from sunlight.Storage procedures are the main factor in determining glove shelf life. Gloves should be kept in their packaging protected from sunlight, artificial light, humidity and stored at temperatures between 10 °C - 30 °C. Real time stability testing for expiration date on this products after accelerated ageing process are declared as 5 years.

UK REP
Kerbl UK Ltd 8 Lands End Way
Oakham, Rutland LE15 6RF UK
www.kerbl.co.uk Email: enquiries@kerbl.co.uk

UK REP
Kerbl UK Ltd 8 Lands End Way
Oakham, Rutland LE15 6RF UK
www.kerbl.co.uk Email: enquiries@kerbl.co.uk

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