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Table of Contents

1. Product Overview 25
1.1 Correct Use 25
1.2 Introduction of main unit 25
2.1 Power supply 25
2.2 Operation environment 25
2.3 Storage and Transportation 25
2.3.1 Storage and transportation environment 25
2.3.2 Storage and Transportation 26
2.4 Further Technical Specifications 26

3. 3 System Configuration 26
3.1 Standard configuration 26

4. System Check 26
4.1 Place the instrument 26
4.2 Check the instrument before and after use 26
4.3 System reset 27
4.4 Charge the instrument 27

5. Getting Started 27
5.1 Preparation before use 28
5.2 Getting Started 28
5.3 Change the APP language 28

6. Functional Operation 28
6.1 Power on/off 28
6.2 Setting 28
6.2.1 Input the barn name 28
6.2.2 Set the image format 29
6.2.3 Set the cine frames 29
6.2.4 Set the probe protect time 29
6.3 Automatic ear tag scanning 29
6.4 New User 29
6.5 Image adjustment and image control 29
6.5.1 Total gain adjustment 29
6.5.2 Near field gain adjustment 29
6.5.3 Far field gain adjustment 29
6.5.4 Dynamic range adjustment 29
6.5.5 Frequency adjustment (Frequency conversion) 29
6.5.6 Depth range selection 29
6.5.7 Image Freeze/Unfreeze 30
6.6 Lock screen 30
6.7 Save image/video 30
6.8 Play back and Measurement 30
6.9 Save check list 30
6.10 Record query 30
6.11 Notes 30

7. Measurement 30
7.1 General measurement 30
7.1.1 Distance measurement (G.A and EDC measurement) 31
7.1.2 Ellipse measurement 31
7.1.3 Punctuation measurement 31
7.2 Measurement of Swine's lean percentage 31

8. System Maintenance 32
8.1 Daily maintenance 32
8.1.1 Cleaning and disinfection 32
8.2 Use and maintenance for the rechargeable battery 32
8.2.1 Use and maintenance for the rechargeable battery 32
8.2.2 The maintenance for the rechargeable battery 32
8.3 Troubleshooting 33

9. Further Notes 33
9.1 Warranty and repair service 33
9.2 Electronic Waste 33
9.3 Declaration of Conformity 33

Spare Parts 88


(EN) Instructions for use
Wireless ultasonic scanner WRLS #221300

1. Product Overview

1.1 Correct Use

The WRLS Ultrasound Device for veterinary use is composed of main unit, ultrasound APP and mobile terminal. Equipped with a 3.5 MHz mechanical sector probe, it is connected to a tablet or smartphone with an ultrasound APP via WiFi and enables ultrasound examination to be performed right on site.

1.2 Introduction of main unit

No	Name	Description
1	Beam focus	To convert electric signal to ultrasonic signal based on principle of converse piezoelectric effect. The ultrasonic signal, after entering the body, is reflected as echo wave and converted to electric signal again. The beam focus is on the probe surface. Supply ultrasonic coupling gel to the ultrasonic transducer when performing ultrasonic diagnosis.
2	Freeze/unfreeze-key	Press the key freeze/unfreeze the image.
3	Power switch	Power on or power off the instrument.
4	Power supply indicator	Turn on the instrument, the indicator light is on, and turn off the instrument, the indicator light is off. The indicator light blinks when the battery is low.
5	Battery charge indicator	When charging, the indicator light is on, and when it is full, the indicator light is off.
6	Charging port	Connect the TYPE C data cable to charge.
7	Ear tag scan position	The position used to scan ear tag.

2. Product Specifications

2.1 Power supply

Power supply input: 100V-240V~, 50Hz/60Hz
Main unit power input: DC5.0V 3.0A
Internal supply voltage: DC3.6V-4.2V

2.2 Operation environment

Ambient temperature: 10°C-40°C
Relative humidity: 30%-75% (without condensation)
Atmospheric pressure: 700hPa-1060hPa

2.3 Storage and Transportation

2.3.1 Storage and transportation environment

Ambient temperature: -20°C-55°C
Relative humidity: 30%-93% (without condensation)
Atmospheric pressure: 700hPa-1060hPa

2.3.2 Storage and Transportation

1. If the instrument is stored over 3 months, take out the instrument from the packing case, connect it to power supply for 4 hours, and then disconnect the power. Place it in the case again following the direction indicated by arrows on the package. Store the case in the warehouse. Do not pile the case. The instrument case should have adequate space from ground, walls and ceiling of the warehouse.
2. Environment requirement Ambient temperature: -20°C-55°C; Relative humidity: 30%-93% (without condensation); Atmospheric pressure: 700hPa-1060hPa. The warehouse should be well ventilated and free of direct sunlight and corrosive gas.
3. Shockproof measures have been taken inside the packing case to allow for transport by air, railway, land and sea. The goods shall not be exposed to poor weather conditions like rain and snow, nor shall the goods be placed upside down, bumped, knocked or over-stacked.

2.4 Further Technical Specifications

1. Main unit size: approx. 231mm (length) ×52mm (width) ×79mm (height)
2. Weight of main unit: approx. 410g
3. Transmission frequency 2.4 GHz, transmission power ≤ 20 dBm (100 mW)

3. 3 System Configuration

3.1 Standard configuration

1	Scanner	1 PC
2	Ultrasound APP software	1 PC
3	Data cable	1 PC
4	Phone wristband	1 PC
5	Lanyard	1 PC

4. System Check

⚠ Warning: The waterproof grade of the instrument is IPX4 (no adverse effect on splashing water in all directions). Splash any liquid on the instrument may damage it. If you accidentally splash liquid on the instrument, please immediately turn off the power and contact your service representative.

⚠ Warning:

1. When the instrument works abnormally, do stop working, turn off the power and check the reason, then contact Kerbl about it.
2. Turn off power after ultrasonic diagnostic operation.

4.1 Place the instrument

Please carefully read through and fully understand the safety cautions before moving and placing the instrument.

1. Unpack the instrument case and check the goods for its completeness.
2. Place the instrument in a stable and leveled position.

4.2 Check the instrument before and after use

1. Before and after ultrasonic diagnosis to check if there are any exceptionally on the beam focus, such as peeling, cracks, bulge, or if the transducer is reliable, cleaned and disinfected.
2. Check the instrument to ensure a normal start-up, and no abnormal voice, strange smell or overheating.
3. Check the transducer for abnormal heat when the instrument is in use. You can touch the probe by hand to feel its temperature.
4. Check the image to ensure trouble-free display (no abnormal noise or flicker).

⚠ Attention: If the overheat transducer is placed on the patient's skin, heat injury may occur.

⚠ Attention: Thoroughly clean the coupling gel on the probe surface each time after ultrasonic operation, or the coupling gel may become hardened on the transducer of the probe, deteriorating quality of image.

4.3 System reset
In case of abnormal screen display or no-working for system operation, try to restart the system by turning on/off the power.

4.4 Charge the instrument
Charge the instrument with the TYPE C charger on the phone/Pad or portable charger. The specifications for the charger are 5.0V, 2.0-3.0A.
Instrument can be used for about 4.5 hours after fully charged; the shortest charging time is 3.5-4 hours, up to about 5 hours; Over-charging or discharging will shorten the battery life. When charging, the battery charge indicator light is on; after fully charged, the indicator light is off.

Note: The device is supplied without its own charger. Standard chargers for smartphones or tablets can be used.

⚠ Attention: When the power is too low, the power indicator light blinks, and the instrument will automatically shut down, and then should be charged.

⚠ Attention: The battery charger should conform to the CCC or the country's phone charger standard requirements.

⚠ Warning:
1. **When charging is complete or not charged, disconnect the charger from the instrument and unplug the charger from the power outlet.**
2. **If the charger plug, shell or power cable has damaged, don't continue to use to avoid electric shock or fire.**
3. **If the charger is wet, soaked or severe moisture, please stops using and contact the Kerbl company service agents.**
4. **Do not use the instrument while charging, may cause a safety risk.**
5. **Use unapproved or incompatible power supply, charger or battery, which may cause a fire, explosion or other hazards.**


5. Getting Started

⚠ Attention: When installing the "WRLS APP", if the installation is blocked, please allow to continue installing it.

⚠ Attention: After the instrument is connected to the phone or pad, if appear "Allow WRLS to access this device's location", "Allow WRLS to access photos, media, and other files on your device", please choose allow; if appear "This WLAN network has no Internet access", please choose "Connect".

⚠ Attention: After the instrument is running for a while, the temperature of the instrument will rise. If the temperature is too high, please do not contact for a long time, otherwise it may lead to low temperature burns, causing skin swelling or pigmentation.

Before using the instrument, you need to install "WRLS APP" on your Pad or Smart Phone. "WRLS APP" is suitable for Android phones or tablets only. WRLS APP supports landscape or portrait display. All the icons on the APP may be different due to the display screen of the phone or tablet; the display interface may be slightly different; the actual products shall govern. The usage of the phone and tablet is similar, this article only take the landscape screen of the phone as an example to introduce.

WLAN

5.1 Preparation before use
1. Obtain “WRLS APP” and install it.

Note: According to the information provided in the APP download method label pasted on the machine, obtain APP software. Remove this label before using it.

Software Download
The required software can be downloaded from the App Store/Google Play or from the website www.kerbl.com/de/produkte/agrarbedarf-9901/veterinaerbedarf-9996/traechtigkeitsdiagnose-10015.

App Download

2. After finishing the installation, set up the phone.
(1) Go to Settings, tap Display, tap Font size to adjust text size on screen, in order to display completely, should select “Normal”.
(2) Go to Settings, tap Permissions to enter permission interface, tap “WRLS”, and enable the “Do Not Disturb” permission so that to adjust the volume of the ear tag scan tone.
3. After finishing the operation, you can use it.

5.2 Getting Started
1. Power on the instrument. Hold down the power switch „⏻” for 3-5 seconds, and release the switch when the power supply indicator light is on;
2. Connect to WLAN. Go to phone Settings, tap WLAN, open WLAN, the phone will search and list the available WLAN;
3. In the available WLAN list, select WRLS network ID, such as “W36666667”, connect to it;
4. After connected successfully, tap the WRLS APP icon to enter the ultrasound interface, prompt “The initialization data ends!”, the wireless signal icon “” on the APP is green, then you can use the APP.

Note:
1. The default WLAN ID of the instrument is the number displayed on the shell of the instrument. For example, “W36666667”.
2. If the connection fails for the first time, exit the APP and re-enter.
3. When the network is disconnected, the wireless signal icon “” on the APP will appear white and blink continuously, and the instrument will automatically turn off
after 3 minutes.


5.3 Change the APP language
1. Go to phone Settings, tap Language & input, tap Language and select the language to use;
2. After the language setting is complete, enter the ultrasound APP again, and the ultrasound interface will display the changed language.

Note: WRLS APP supports Chinese, English, German, French, Spanish and Polish. Other languages are not supported.

6. Functional Operation

6.1 Power on/off
Power on: Hold down the power switch ⏻ for 3-5 seconds, and release the switch when the power supply indicator light is on;
Power off: Hold down the power switch ⏻ and release the switch when the power supply indicator light is off.

6.2 Setting
In the frozen state, click Setting key to display Setting dialog box, as shown below. In the Setting dialog box, you can set barn name, image format, cine frames, probe protect time also, can view the video format, software full version.

Note: The battery level and time displayed on the APP are synchronized with the phone/pad.

6.2.1 Input the barn name
1. In the Setting dialog box, tap the barn name input box, the onscreen keyboard will automatically open;
2. After the input is complete, close the keyboard;
3. Click Cancel, exit the setting dialog box;
4. Click Save, the ranch name is displayed on the screen.

6.2.2 Set the image format
1. In the Setting dialog box, select the image format bmp, jpg or png;
2. Click Cancel, exit the setting dialog box;
3. Click Save, confirm the setting and exit the setting dialog box.


6.2.3 Set the cine frames
1. In the Setting dialog box, select the cine frames 256, 512 or 1024;
2. Click Cancel, exit the setting dialog box;
3. Click Save, confirm the setting and exit the setting dialog box.

Note: The cine frames can be selected depending on the mobile terminal configuration. The higher the configuration, the larger the number of selected cine frames.

6.2.4 Set the probe protect time
1. In the Setting dialog box, tap the probe protect time input box and input the probe protection time, the probe protection time can be set within "0 ~ 99" „ minutes;
2. Click Cancel, exit the setting dialog box;
3. Click Save, confirm the setting and exit the setting dialog box.
Note: In real-time state, when the set probe protective time is reached the instrument automatically enters the frozen state. Press freeze/unfreeze key to return to the working state.

6.3 Automatic ear tag scanning
After the instrument is powered on, hold down the freeze/unfreeze key on the instrument for 4 to 6S and keep it pressed, and move the ear tag scan position close to the ear tag. When you hear the prompt tone, the user ID is automatically recorded and displayed on the screen and then release the freeze/unfreeze key.

6.4 New User
1. In the frozen state, click  key to display User dialog box, as shown below;
2. Enter the ID. You can enter the ID in two ways:
 - The ID entered by the scanning ear tag is automatically displayed in the input box (see 6.3).
 - To enter an ID manually:
 (1) Tap "New", and the ID will be automatically generated based on the system date;
 (2) Tap the ID input box and the on-screen keyboard will open automatically, enter the user ID (up to 15 digits);
3. Enter user weight (1-699) kg;
4. Choose male or female;
5. Click Cancel, exit the User dialog box;
6. Click Save, save the user information, and ID, weight are displayed on the screen.

Note: The user information must be all entered, otherwise you cannot click Save.

6.5 Image adjustment and image control

6.5.1 Total gain adjustment
In real-time state, click the arrow on either side of the "Gain" key to adjust the total gain of the image. Click the right arrow to increase the total gain; click the left arrow to decrease the total gain. Control the total gain to get a satisfactory image. The total gain value is displayed on the left side of the Gain key.

6.5.2 Near field gain adjustment
In real-time state, click the arrow on either side of the "Near" key to adjust the near field gain of the image. Click the right arrow to increase the near field gain; click the left arrow to decrease the near field gain, to control the near field gain. The near field gain value is displayed on the left side of the Near key.

6.5.3 Far field gain adjustment
In real-time state, click the arrow on either side of the "Far" key to adjust the far field gain of the image. Click the right arrow to increase the far field gain; click the left arrow to decrease the far field gain, to control the far field gain. The far field gain value is displayed on the left side of the Far key.

6.5.4 Dynamic range adjustment
In real-time state, click the arrow on either side of the "Dync" key to adjust the dynamic range of the image. Click the right arrow to increase the value of dynamic range; click the left arrow to decrease the value of dynamic range, to control the dynamic range of the image. The dynamic range value is displayed on the left side of the Dync key.

6.5.5 Frequency adjustment (Frequency conversion)
In real-time state, click the arrow on either side of the "Freq" key to convert the frequency. Frequency conversion is 2.5M/3.5M. The frequency value is displayed on the left side of the Freq key.

6.5.6 Depth range selection
In real-time state, click the arrow on either side of the "Depth" key to get images in different depth. There are 7 kinds of depth range. The depth value is displayed on the left side of the Depth key, as well as at the "D:" in the lower left corner of the image area.

6.5.7 Image Freeze/Unfreeze

There are two ways to implement:

1. In real-time state, click key to freeze the image; in frozen state, click key to unfreeze the image.
2.Press key on the instrument to freeze or unfreeze the image.

6.6 Lock screen

Click the icon in the upper left corner of the screen to appear the lock screen icon "", the screen is locked and cannot be turned over automatically. Click the icon again to exit the lock screen.

6.7 Save image/video

1. Freeze the current image, click "Save" key to display Save dialog box;
2. Tap the text input box, the cursor is in the input box, the onscreen keyboard will automatically open and enter the file name;
3. Click Cancel, exit the Save dialog box;
4. Click Save Image or Save Video, the screen appears "Save Success" prompt.

Note: The saved images/videos are stored in the folder of the phone or pad by default. Follow following path:

6.8 Play back and Measurement

Play back

In real-time state, the system is always saving the scanning image. The playback images are for a period time images before freezing.
Freeze image, click key to enter the automatic playback status, and then click key to return to manual playback status, click keys to choose the image needed playback. Press key unfreeze and exit playback status.

Playback measurement

In manual playback status, choose measured image and perform measurement, the detail measurement method refers to corresponding chapters.
Press key again to unfreeze and exit playback status.

Note: If the images appear abnormal, that is without enough storage time, and the images have not been stored full.

6.9 Save check list

1. Freeze the desired image, click "List" key to display Check List dialog box;
2. Select the desired status Pregnancy/Unpregnancy/Suspected/Disease;
3. Click Cancel, exit the Check List dialog box;
4. Click Save, the system saves the user status and displays it on the screen.

6.10 Record query

1. In the frozen state, click "Record" key to display Record screen;
2. In the "Search" bar, set the search criteria: ID, status, and check date;
3. Click Search, the eligible records will be displayed in the check list and file list below;
4. Click Page Up, Page Down to turn the page to view the records; you can directly click the images or videos in the file list to view;
5. Click Clear to clear all records in the search bar, check list, and file list;
6. Click Exit to Exit the record screen.

Note:

1. The check list records the status of each check. The list is as follows: ID, sex, status, weight, and check date.
2. The file list records the saved images and videos during each check.

6.11 Notes

1. In the frozen state, click "Note" key, this key is highlighted in green;
2. Tap the image area where need to be annotated, the cursor "| " appears on the screen and the onscreen keyboard will automatically open;
3. After entering the annotations, move the cursor position to make multiple notes;
4. If you want to modify the contents during the note, move the cursor "| " to the character need to be modified, delete the character or re-enter it;
5. Click "Note" key again to complete the note entry;
6. Click key to clear the notes;
7. Press key to unfreeze and clear the screen.

7. Measurement

7.1 General measurement

In the frozen state, click "Meas" key to pop up measurement menu, showing three measurement methods Distance, Ellipse or Punctuation. Click the different measurement methods to enter the corresponding measurement interface, as shown below:

Discance
Ellipse
Punctuation

⚠ Attention: The accuracy of software measurement: distance measurement ≤ 0.1mm; area measurement ≤ 0.01cm²; volume measurement ≤ 0.01cm3. Due to differences in images obtained by each user in different times, the actual object for the accuracy of the measurement may be greater than the above-mentioned values.

⚠ Attention: When measuring, the positioning icon appears on the screen. Hold down any four corners at the rectangular frame of the positioning icon to slide, you can move the positioning icon to anywhere. Hold down the direction keys of the positioning icon to fine-tune the position of measurement mark, tap the middle Confirm key to end the current measurement.

⚠ Attention: When the display screen is small, the measured results displayed on the right side of the screen may not be displayed completely. In this case, you can slide up or down, or hide virtual keys on the phone to view.

7.1.1 Distance measurement (G.A and EDC measurement)

Function: measure the gestational sac size of swine and calculate the gestational age (G.A) and estimated date of confinement (EDC).

Operation steps:

1. Freeze desired image, click "Distance" item in the measurement menu to start distance measurement, and the positioning icon appears on the screen at the same time;
2. Tap the image with your finger, appear white mark "✕",move your finger to appear measuring line and green mark "✕", move the green mark to locate the starting point of the measurement;
3. Tap white mark "✕" of the other endpoint, then the white mark turns green, move the green mark to the endpoint of the measurement, the measured distance value real-time shown on the right side of the screen;
4. Click the "✕"mark to exchange the starting and end point; Hold down the direction keys of the positioning icon to fine-tune the position of starting and end point;
5. After completing the measurement, tap the Confirm key mark to exchange the starting and end point; Hold down the direction keys of the positioning icon to fine-
 tune the position of starting and end point;
6. Click key to clear measurement marks and the results.

7.1.2 Ellipse measurement

Function: measuring the volume of the subject. 2-Axis method only needs to measure the vertical section of the subject.

Formula of 2-axis method: $(\pi/6) * A * B^2$

Where, A is the major axis of the ellipse; B is the minor axis of the ellipse. Operation steps:

1. Click "Ellipse" item in the measurement menu to start ellipse measurement, and the positioning icon appears on the screen at the same time;
2. Move your finger to display ellipse and green mark „✕" move the green mark to locate the starting point of the measurement;
3. Tap white mark "✕" " " of the other endpoint, then the white mark turns green, move the green mark to the endpoint of the measurement, the measured circumference, area and volume values real-time shown on the right side of the screen;
4. Tap the measurement mark "✕" on the elliptical arc, adjust the elliptical minor axis;
5. Repeatedly tap the measurement mark "✕", you can exchange the location of the measurement points; hold down the direction keys of the positioning icon to fine-tune the position of measurement mark;
6. After completing the ellipse measurement, tap the Confirm key on the positioning icon, the measured circumference C, area A and volume V values are displayed on the right side of the screen, and the measurement is finished;
7. Click key to clear measurement marks and the results.

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1. Click "Punctuation" item in the measurement menu to start punctuation measurement, and the positioning icon appears on the screen at the same time;Tap the starting point of the measurement with your finger, appear white mark "✕", the location of the white mark as the first point of the measurement; Along the boundary of the measurement range, continue to tap with your finger to mark the required points;
2. When completed the desired punctuation, tap the Confirm key on the positioning icon, appear the trace on the screen (trace is automatically drawn along the order of marking point until close the beginning point and end point), the measured circumference and area values automatically shown on the right side of the screen; If continuously marks the thirty-two point, it directly appears trace on the screen and the measured circumference C and area A values are automatically displayed on the right part of the screen;
3. Click key to clear measurement marks and the results.

7.2 Measurement of Swine's lean percentage

% Lean calculation uses NSIF formula.

1. In the frozen state, click "Lean" key to start lean percentage measurement, a middle vertical line and the positioning icon appears on the screen at the same time;
2. Measure fat thickness. Slide the green mark "✕" with your finger to the endpoint, now the mark can only move vertically, the fat thickness value is displaying on the right side of the screen in real-time, tap the Confirm key , the fat thickness value is displayed on the right side of the screen;
3. Measure loin. The measurement method refers to the "Distance" measurement of general measurement;
4. After completed the fat thickness and loin measurement, the result of %Lean is automatically displayed on the right of the screen;
5. Click key to clear measurement marks and the results.

Note:

Before measuring the lean percentage of the swine, you must enter the swine's weight (range 1 ~ 699kg) otherwise, the default weight of swine is 100kg.

8. System Maintenance

The system maintenance should be performed by the user. Users shall be in full charge of maintenance and operation of the system after purchasing the product. Please check the status of the machine regularly and perform a preventive check for the system.

⚠ Warning: The following maintenance must not be performed when inspecting.

8.1 Daily maintenance

Daily maintenance should be performed by the user.

8.1.1 Cleaning and disinfection

⚠ Warning: Turn off the instrument before cleaning. Cleaning in the power-on may cause the instrument to be damaged.

⚠ Warning: The waterproof grade of the instrument is IPX4 (no adverse effect on splashing water in all directions). Do not splash any water or liquor into the instrument when cleaning or maintaining otherwise, it will cause malfunction.

⚠ Attention:

1. To prevent possible infection, it is advisable to wear sterilized gloves when cleaning, disinfecting the instrument.

2. In the process of cleaning and disinfection, avoid probe overheat (exceeding 60°C) as it may be deformed or damaged under excessive heat.

1. Clean the instrument
 (1) Must wear sterilized gloves to prevent possible infection.
 (2) Wipe the probe with water or soapy water to remove all contaminants. Do not use brushes as it may damage the probe.
 (3) After finishing the wiping, use a sterilized cloth or gauze to wipe the water on the probe. Do not dry the probe by heating it.

2. Disinfect the instrument
 (1) Routine disinfection must wear sterilized gloves to prevent possible infection;
 (2) Clean the probe before routine disinfection, and then use 75% alcohol solution to wipe the probe twice;
 (3) Use sterile water to rinse the remaining chemicals on the probe;
 (4) Use sterilized cloth or gauze to wipe the water on the surface of probe. Do not dry the probe by heating it.

⚠ Attention:

1. It is a normal phenomenon that the beam focus may discolor and the text on the instrument label may fade away.

2. The regular disinfection times should be minimized as it may lead to degrade of the probe safety and performance.

3. Check instrument after cleaning and disinfection
 (1) Check the instrument enclosure to ensure they are free of abnormity (such as scuffing, cracks or drop-off);
 (2) Check the transducer of instruments to ensure that there are no any abnormity, such as scuffing, cracks, peeling, bulge and drop-off.

8.2 Use and maintenance for the rechargeable battery

8.2.1 Use and maintenance for the rechargeable battery

1. Only use battery provided by Kerbl Company.
2. Do not remove the built-in rechargeable battery to avoid damage to the battery or the instrument. The battery should be replaced only by Kerbl company authorized service personnel.

8.2.2 The maintenance for the rechargeable battery

1. Battery is consumable; the battery cycle-life is based on the times of charge and discharge as unit. When the use time reduced significantly compared with normal conditions, the battery should be promptly replaced.
2. The excess high or low temperature will affect the charging and discharging performance, and shorten the battery life and capacity.

⚠ Attention: The battery charger should conform to the CCC or the country's phone charger standard requirements.

⚠ Attention: Don't throw away the exhausted battery anywhere; especially throw it in the fire. Please deal with it according to local statutes. Use pollution degree II to deal with.

Attention:

1. **Do not throw the battery into water, which will lead to the battery leakage, explosion or fire;**
2. **Do not use or store the battery near the heat source, such as fire or heater, which will lead to the battery leakage, explosion or fire;**
3. **Do not connect the anode and cathode reversely, which will lead to battery leakage, explosion or fire;**
4. **Do not heat up or throw the battery into fire, which will lead to leakage, explosion or fire;**
5. **Do not connect the anode and cathode with any metal or conductor; do not transport or store the battery together with necklaces, hairpins or other metal objects, which will lead to leakage, explosion or fire;**
6. **Do not hammer blow, throw or mechanically shake the battery, which will lead to leakage, explosion or fire;**
7. **Do not insert the battery with nail or other sharp objects; do not hammer or trample the battery, which will lead to leakage, explosion or fire;**
8. **Do not weld the battery terminal directly, which will lead to leakage, explosion or fire;**
9. **Do not disassemble the battery in any way, which will lead to leakage, explosion or fire;**
10. **Do not charge the battery near the heat source or extra-hot environment, which will lead to leakage, explosion or fire;**
11. **Do not put the battery into the microwave oven or pressure vessel, which will lead to leakage, explosion or fire;**
12. **Do not mix the battery together with one-off battery (such as dry battery), or different capability or different model or different brand battery, which will lead to leakage, explosion or fire;**
13. **Do not use the abnormal battery with smell or abnormal heat or distortion or turn colors or abnormal phenomena, which will lead to leakage, explosion or fire;**
14. **Do stop the charge and unplug from the charger at once if any abnormal phenomenon happens to the battery, such as smell or abnormal heat or distortion or turn colors. Otherwise, each of above will lead to leakage, explosion or fire;**
15. **Do remove the battery from the nearby fire if any leakage or smell happens, which will lead to the leakage, explosion or fire;**
16. **If any leakage splash into eye, do not wipe the eye, instead of washing it and get help from the doctor as soon as possible. Otherwise, the eye will be injured;**
17. **Do not use the battery in an extremely hot environment, such as hot sunshine or in the car when it is too hot, because these will catch fire, even worsen its performance and shorten its life;**
18. **If use the battery beyond the listed environment on the manual, it will worsen its performance or shorten its life, even lead to extreme heat or explosion or fire.**

8.3 Troubleshooting

To ensure normal operation, users are recommended to prepare a proper maintenance plan.
If the following problems occur after turning on the instrument, try to make corrections following the methods in the table. If the problem remains unsolved, please contact Kerbl company for support.

SN	Troubles	Correction
1	Power supply indicator does not light when powering on the machine.	1. Check power supply. 2. Check whether the instrument is powered on.
2	Character and gray scale are displayed, but no ultrasonic image on the ultrasound APP.	1. Check whether the instrument is connected to the WLAN. 2. Turn on the instrument, connect to WLAN, wait 5-6 seconds and then open the ultrasound APP.
3	The instrument is not connected to WLAN.	1. Check whether the mobile terminal device is opened the WLAN. 2. Check that the mobile terminal device is properly connected to the WLAN.

9. Further Notes

9.1 Warranty and repair service

Albert Kerbl GmbH grants you a 2-year guarantee on this article in accordance with our general guarantee conditions.
These general warranty conditions can be found at www.kerbl.com/guarantee. If the device is used in a manner that does not comply with the intended use, then this voids the manufacturer's warranty and liabilities. All repairs are all performed by professionals recognized by Kerbl; all replacement parts and accessories, consumables involved repairs are Kerbl company.

9.2 Electronic Waste

Do not dispose of batteries in domestic waste. Dispose of the battery properly, separately from the electrical device and in accordance with locally applicable regulations. Municipal collection points, as well as return points operated by retailers, can be used to return such batteries. The separate collection and reprocessing of old batteries is intended to avoid negative impacts on the environment and human health.

9.3 Declaration of Conformity

CE

Albert Kerbl GmbH hereby declares that the device described in the instructions complies with Directive 2014/53/EU. The full text of the EU Declaration of Conformity is available at the following web address: www.kerbl.com/doc

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Kerbl UK Ltd 8 Lands End Way