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2. Installation 30
2.1 Ambient conditions 30
2.2 Installation and installation of the probe 30
2.3 Power supply 30
2.4 Battery indicator 31

3. Operation 31
3.1 Display 31
3.2 Control panel 31

4. Operation 32
4.1 Switching on 32
4.2 Examination and Diagnosis 32
4.3 Settingimage parameters 33
4.3.1 Setting the frequency 33
4.3.2 Adjusting the gain 33
4.3.3 Adjustingbrightness and contrast 33
4.4 Input menu 33
4.5 Image Management 34
4.5.1 Saving an image 34
4.5.2 Viewing Images 34
4.5.3 Image Smoothing 35
4.5.4 Colour Mode 35
4.5.5 Circumference , area and volume measurement 35
4.5.6 Histogram statistics function 36
4.6 Distance measurement 36
4.7 Volume measurement 37
4.8 Heartbeat measurement 37
4.8.1 Heart rate measurements using the spacer method 38
4.8.2 Heart rate measurements using the circumference method 38
4.9 Estimating the gestation period in pigs 38
4.10 Other 38
4.10.1 Printing 38
4.10.2 Transferring images to a PC 38
4.10.3 Switching off the device 38

5. System maintenance 38
5.1 Daily maintenance 39
5.2 Cleaning and Disinfection 39
5.3 Use for and replacement of the rechargeable battery 39
5.4 Maintenance of the rechargeable battery 39

6. General instructions 40
6.1 Warranty and Repairs 40
6.2 Electronic Scrap 40
6.3 CE Declaration of Conformity 40

Spare part 119

(EN) Instructions for use
Ultrasound scanner including PRO probe #221301

1. Product overview

1.1 Intended Use

The PRO ultrasound scanner is a portable ultrasound device for veterinary applications. It consists of a main unit and a wired 3.5 MHz sector probe and is intended for ultrasound examinations of pigs, sheep, goats and alpacas.

1.2 Device Overview


Fig. 1: Front view


Fig. 1: Side view

1.3 Technical Data

probe		3.5 MHz mechanical sector
Displayed depth (mm)		60–190
Maximum scanning depth (mm)		≥140
Resolution (mm)	Lateral	≤4 (depth ≤80) ≤5 (8≤depth ≤130)
	Axial	≤2 (depth ≤80)
Blind zone (mm)		≤8
Geometric	Horizontal	≤20
Position accuracy (%)	Vertical	≤10
Monitor size (inches)		5.6” TFT-LCD
Display modes		B, B+B, B+M, M
Greyscale		256 levels
zoom		x1.0, x1.2, x1.5, x2.0
Image storage		128 frames
Measurement functions		Distance, circumference, area, heart rate, stage of pregnancy
battery capacity		3350mAh
Device power		13W in use / 25W whilst charging
Charger power consumption		45W
Data display		Version date, age
Colour display		Grey, red, yellow, blue

1.4 Basic Principle
B-mode ultrasound provides visibility into internal structures by transmitting ultrasound waves into the tissue and analysing the returning echoes. This works as follows: Different types of tissue have different densities and therefore different levels of resistance to ultrasound, i.e. different acoustic impedances (a combination of sound density and sound velocity). When the piezoelectric chip (transducer) receives a specific electrical impulse, it generates ultrasound at a specific frequency. When these ultrasound waves penetrate the tissue, the various organs produce different echoes. The differently reflected sound waves are received again by the transducer, which converts them into electrical impulses, reinforces them, organises them and digitises them. A standard video signal is generated and cross-sectional images of organs are displayed on the screen.

1.5 Standard configuration

221301-1	PRO Main Unit
221301-2	3.5 MHz mechanical sector probe
221301-3	Li-ion battery
221301-4	Neck strap with splash-proof cover
221301-5	Bottle holder for ultrasound gel
221301-6	Carry case
221301-7	Charger (EU) for PRO
221301-8	UK plug with cable for PRO

2. Installation

2.1 Ambient conditions
1. Ambient temperature range: +10 °C to +40 °C
2. Ambient humidity: 30 % to 75 %
3. Atmospheric pressure range: 70 kPa to 106 kPa


Avoid severe vibrations during operation and keep the device away from equipment generating strong electric fields, strong magnetic fields or high voltage. Do not expose the display to direct sunlight. Ensure the device is well ventilated and provide it with protection against moisture and dust.

2.2 Installation and installation of the probe
There is only one probe socket to which all supported probes can be connected.

⚠ Warning: Under no circumstancesshould you disconnect or reconnect the probe plug whilst the device is switched on if the probe or the main unit is damaged. Once the probe is connected to the main unit, do not simply disconnect or reconnect it, even if there are connection problems.

⚠ Important: The transducer is delicate; provide protection from falling and breaking. Handle the device with care.

2.3 Power supply
1. Insert the battery correctly into the main unit: Remove the battery compartment cover by pushing the squeegee towards the front. Insert the battery.

2. Press the power switch on the main unit to switch it on; the power indicator will then light up.
3. The device is now ready for use.

⚠ Note: If the main unit runs out of power, a beep will sound, indicating that the battery is empty and needs recharging.

⚠ Warning: You must not use any power pack other than the one supplied.


2.4 Battery indicator

There is a battery indicator in the top right-hand corner of the screen. When the main unit runs out of power, a beep will sound, indicating that the battery is empty and needs recharging.


1. Plug the DC14V/3A output of the adapter into the DC14V charging socket on the side.
2. Plug the 'Input AC220V, 50Hz' power cable into a socket.
3. Regardless of whether the main unit is switched on or off, if the 'Low Voltage' indicator is flashing, the charger is charging the battery. When the 'Low Voltage' indicator goes out, the battery is fully charged.

3. Operation

3.1 Display

3.2 Control panel

Freeze image

Press this button to switch between FREEZE and live view

Tip: 'FREEZE' appears in the bottom right-hand corner of the screen, indicating that the image is frozen.

Function buttons:

M/O M-mode display button

Press this button in B/M mode whilst in real-time or freeze mode; the display will now show an image in M mode.

Press this button outside B/M mode to open the image editing menu. Press CLEAR to exit this function.

1/ZOOM Set scan depth

If you are in B-mode during a measurement, you can use this button to adjust the scan depth (adjustable range from 60 to 190 mm, in 14 steps).

A typical scan depth for pregnancy testing in pigs is 160 mm

+/2 Setting measurement points

If you are in the FREEZE screen, you can use this button to set a '+' marking, which you can move using the direction buttons. (For detailed instructions, see 'Instructions for measuring distance, circumference and area')

▲3 ~ ▼6 Arrow keys

If you are in a measurement mode, press these buttons to move the cursor.

If you are in the B-mode LIVE image, you can vary the scan depth by pressing the ▲3 and ▼6 buttons. The current scan depth is displayed below the screen.

In the input menu under "name" and "pasture", press ▲3 and ▼6 to display the desired settings.

Press 5◀ or ▶4 to adjust the parameters on the right-hand side. The parameter currently available for adjustment will be highlighted. Press ▲3 or ▼6 to adjust it.

B/7 B-Mode Display Button
Press this button whilst in LIVE or FREEZE view to enter B-mode (Single B-mode is the default setting).

BB/8 BB Mode Display Button
Press this button in the LIVE or FREEZE view to enter BB mode. A LIVE view and a FREEZE view are displayed one-sidedly. Press repeatedly to switch between the FREEZE view and the LIVE view. If you press , two FREEZE views are displayed.

BM/9 B/M mode display button
Press this button twice in the LIVE or FREEZE view to enter B and M mode (B/M mode). The screen displays B/M mode. The left-hand side shows the B-mode LIVE view, whilst the right-hand side shows the M-mode LIVE view.
Press repeatedly to switch between B/M and M modes in the LIVE image.
In the B-mode image, there is a vertical line (pattern line) showing a number of points with regular spacers. Press 5◀ to move the pattern line to the left. Press ▶4 to move the pattern line to the right.

M/O M Mode button
If you are in B/M mode, you can use this button to switch to the M-mode display status in the LIVE or FREEZE image. The display shows an M-mode image.
If you are not in B/M mode, you can use this button to display the image processing menu. Press the CLEAR button to exit the function.


Fig. B, BB, B/M and M modes (from left to right)

LIST Calving supplies menu
Press this button in the FREEZE screen in B or BB mode to display the list menu. Press the number keys and carry out the relevant operations. For details on operation, refer to the instructions below. Press CLEAR to exit the list menu.

DISTANCE Operation key for measurements
Use ' +/2 ', ' M/O ' and the directional keys to perform distance, circumference, area, volume and histogram measurements. For details on operation, refer to the instructions below.

OK Menu button Comment
Press this button to display the input menu. Press the number keys to access the sub-menu. Press to exit the input menu.

CLEAR Reset screen settings
Pressing this button clears the measurement mark, adds a comment to results or exits the relevant menu.

RESET Reset button
If the device stops responding due to an unexpected situation or malfunction (i.e. it no longer responds when any button is pressed), press this button. The device will be reset.

Note: The buttons ' M/O ', ' ⊕/1 ', ' +/2 ', ' B/7 ', ' BB/8 ' and ' BM/9 ' are multifunction buttons, meaning they can perform different functions depending on the application. On the one hand, they are used to select the various modes; on the other hand, they can be used to select letters in the comment function.

4. Operation

4.1 Switching on
Switch on the power switch; the power indicator on the control panel will light up and the start-up screen will appear. Press any button (except and CLEAR to display the scan status. Adjust the display brightness, contrast and overall gain control to achieve a visually satisfactory display.

Note: The image appearance on an LCD display can vary significantly depending on the viewing angle. View the screen straight from the front as possible to best assess brightness, contrast and image detail.

⚠ Warning: Do notcover the ventilation openings. Otherwise, the appliance may overheat and be damaged.

4.2 Examination and Diagnosis
Apply a sufficient quantity of medical ultrasound gel to the probe and place it on the area of the animal's body to be examined. The ultrasound image will then be displayed on the screen. Move the probe slowly until the desired examination area is optimally visualised. If necessary, adjust the overall gain, near-field and far-field gain, and the focus to achieve the best possible image quality.

⚠ note:
1. Apply only moderate pressure when pressing the probe against the examination area. Excessive pressure may damage the probe or cause the animal unnecessary pain.
2. Use a probe and frequency suitable for the examination in question.

4.3 Settingimage parameters

Adjusting the image parameters involves adjusting the probe frequency, focus, frame, image smoothing, gamma calibration, gain, brightness and contrast.
In the LIVE view, press ' 5◀ ' or ' ▶4 '; the currently selected parameter is highlighted with a coloured background. Parameters can be set via ' ▲3 ' and ' ▼6 ' and are displayed in the right top corner of the screen.

4.3.1 Setting the frequency

In the LIVE view, press 5◀ or ▶4 and select the 'Frequency' parameter. Adjust the frequency using ▲3 and ▼6 ; the available options are 2.5 MHz, 3.5 MHz and 5.0 MHz.

4.3.2 Adjusting the gain

In the LIVE image, press 5◀ or ▶4 , select the "Overall Gain" parameter and adjust the value using ▲3 and ▼6 . The adjustment range is from 0 to 60.

4.3.3 Adjustingbrightness and contrast

Press 5◀ or ▶4 , highlight ☼ and ◐ in the top right-hand corner of the screen, and then adjust the brightness and contrast using ▲3 and ▼6 .

4.4 Input menu

Press OK on the FREEZE screen; the display will show the input menu as shown in fig. 1. This menu gives you access to functions such as Name, Age, Gender, Comment, Time, Mode, Language, Delete function and beep settings.

V2.03
0. NAME
1. AGE
2. GENDER
3. COMMENT
4. TIME
5. HOSPITAL
6. DELETE
7. DEFAULT SET
8. BEEP-ON

Press M/0 to select '0. NAME' and enter a name or designation. The input screen is as follows:

PLEASE ENTER NAME: •

0-A 1-B 2-C 7-D 8-E 9-F

Below the input window are 26 characters and the space bar. Press ▲3 and ▼6 to display other characters. Pressing the number buttons displays the numbers; up to 15 characters can be entered. Press CLEAR to delete incorrect characters. Once you have finished entering the text, press OK to exit.

Press 1/ZOOM, select '1.AGE' and enter the patient's age. Up to 3 numbers can be entered. The input screen is as follows:

PLEASE ENTER AGE: •

Once you have finished entering the information, press OK to exit.

Press +/2 , select "2.SEX" and enter the patient's gender. The input screen is as follows:

PLEASE ENTER SEX: 1. M 2. F

Press 1/zoom, select "1.M". Press +/2 , select "2.F".

Press ▲3, select "3.COMMENT"; the screen will display a white cursor, followed by 10 comment characters. Press LIST and DISTANCE to display other characters. Press CLEAR to delete characters. Once you have finished entering data, press OK to exit.

Press ▶4 , select "4.TIME" and change the system time and date. The input screen is as follows:

YY-MM-DD HH-MM-SS

Example: If the date and time are 2010-10-9 and 15:30:28, the entry should be:

YY-MM-DD 101009 HH-MM-SS 153028

Once you have finished entering the data, press OK to exit. If you wish to cancel the entry, press OK.
Press ' 5◀ ', select '5. HOSP' and enter the facility. The input screen is as follows:

PLEASE ENTER HOSP:
•

Press ▼6 , select "6.ERASE" and erase the total memory. The input screen appears as follows:

ERASE ALL STORAGE
1. YES 2. NO

Press 1/ZOOM to delete all images. During deletion, the screen displays the status 'ERASING...' in the top-left corner; no inputs or actions are possible whilst this is happening. When the status disappears from the screen, this indicates that the memory has been erased.
Press +/2 to stop the erasure.

⚠ Important: Whilst 'ERASING…' is displayed, the erasure process is not yet complete. No input is possible during this time. Please wait until the process is complete and do not switch off the device to avoid damage or data corruption.

Press B/7 to select "7. DEFAULTSET" and restore the factory settings. All settings will be reset to their factory defaults. Please take care!
Press BB/8 , select "8.BEEP-ON" and switch the beep sound ON or off.
NOTE: The device should not allow the entry of numbers whilst entering NAME(0), COMMENT(3) and Operation (5). Press OK to exit the comment table if required.

4.5 Image Management
In the FREEZE screen, press M/O to display the image management menu. Functions for saving, deleting, smoothing and changing colours are available. Furthermore, area and volume calculations can be performed, and histograms can be displayed.

V2.03
0. SAVE
1. SVLOAD
2. IMAGE-SMOTH
3. COLOUR
4. AREA-VOLUME
5. HISTOGRAM

4.5.1 Saving an image
This device can permanently store up to 128 image files.
If you are satisfied with your FREEZE image, press M/O to display the image management menu. Then press 1/ZOOM to save the current FREEZE image. At the same time, the screen displays the status with the current image number, e.g. "SAVING……05". Once saving is complete, the status message and image number will disappear. Press ⊠ to return to the LIVE image.
Up to 128 images can be saved and automatically numbered based on the time they were saved. If images 01–20 already exist, the next image will be numbered 21. When the memory is full (with a max. of 128 images), a note appears on the screen:

STORAGE IS FULL ERASE NO.01?
1. YES 2. NO

Do you wish to overwrite image no. 01 with the current FREEZE image? Press 1/zoom for 'yes' or press for 'no'.


4.5.2 Viewing Images
In the LIVE or FREEZE image, press M/O to open the image management menu. Then press 1/zoom to select "SVLOAD". The screen will display the following:

PLEASE ENTER STORAGE NO.:
•

Enter the image number in the format '01' and then press OK (if you enter the wrong number, press CLEAR to delete the characters single-handedly, and then re-enter it). Check the first image; the left-bottom screen displays '01/128'. '01' represents the current image number; '128' indicates the storage capacity in terms of the number of images. Press ▲3 or ▼6 to view the images single-handedly.
Press ⊠ to return to the LIVE image. Repeat the above procedure to view further images.

4.5.3 Image Smoothing

In the LIVE view, press M/O to display the image management menu. Then press +/2 to select image smoothing. Repeat this process to change the level of image smoothing. The change is shown in the top right-hand corner of the screen. The setting range comprises the levels 1M0, 1M1, 1M2, 1M3.

Standard image

Smoothed image

4.5.4 Colour Mode

In the LIVE image, press ' M/O ' to access the image settings. Then press 'Colour Mode' on the right-hand side of the screen. Press ' 5◀ ' or ' ▶4 ' to select the parameters, and press ' ▲3 ' and ' ▼6 ' to adjust the parameters.
IC0, IC1, IC2, IC3 for false colours in grey, red, yellow and blue.
CC0, CC1 for character colour in grey and yellow.
BC0, BC1 for background colour in grey and blue.

4.5.5 Circumference , area and volume measurement

In the FREEZE screen, press M/O to access the image settings. Press ▶4 , select 'Area and Volume', and the screen will look as follows:

ENTER:
1. FREEHAND
2. ELLIPSE

Freehand measurement:

Operation
1. Press 1/zoom to select the freehand measurement method. A cursor will appear on the screen. Press the directional buttons to move the cursor to the starting point of the area to be measured.
2. Press DISTANCE and use the directional buttons to move the cursor along the edge to the end point of the area being measured.
3. Press DISTANCE again to complete the circumference and area measurement.

When measuring circumference and area, press ' M/O ' and ' ▶4 '. Alternatively, press ' +/2 ' directly and repeat steps 2 to 3. Up to two sets of data can be measured. The result is displayed on the right-hand screen.



Fig. 1: Measuring circumference and area (freehand measurement method)

C1 and A1 are the first set of circumference and area data.
C2 and A2 is the second set of data for circumference and area.
C1/C2 is the ratio of two perimeter values.
A1/A2 is the ratio of the two area values.

Ellipse measurement:
1. Press and select the ellipse method. The screen will display the ellipse pattern frame.
2. Move the sample frame to the area to be measured.

3. Press DISTANCE and adjust the size using the directional buttons.
4. Press DISTANCE again and use the directional buttons to change the size of the pattern frame.
5. Press DISTANCE and repeat steps 2–4 to readjust the pattern frame.

If you wish to proceed with measuring the circumference and areas, press ' M/O ' and ' ▶4 '. Alternatively, press ' +/2 ' directly and repeat steps 2 to 3. Up to 2 data sets can be measured. The result is displayed on the right-hand screen.

4.5.6 Histogram statistics function

In the FREEZE screen, press M/O to display the settings.
Press 5◀ to select '5.HISTOGRAM'. A selection frame will appear in the image area. Use the directional buttons to select any section of the image. Press DISTANCE to obtain the complete statistics. The result is displayed on the right-hand screen as follows:


Fig. 1: Histogram statistics

The X-axis shows the grey levels. The Y-axis shows the number
PT is the total number of pixels.
Gm is the number of grey levels when the value on the Y-axis is at its maximum.
Pm is the number of Gm pixels
From this image, the total number of pixels (PT) is 10,000; at Gm = 52, the number of pixels is large, at 327.

During operation, press 1/ZOOM or +/2 to zoom in or out on the screen. Press CLEAR to exit the menu.

4.6 Distance measurement

Button settings
1. In freeze mode, press +/2 ; the screen will display the cursor.
2. Use the directional buttons to move the cursor to the starting point to be measured.
3. Press DISTANCE to set the starting point.
4. When you press the directional buttons, another cursor appears. Move it to the end point, press 1/ZOOM to complete the measurement. (Tip: Press DISTANCE in succession to switch between the start and end points)

If you wish to continue measuring distances, repeat steps 1–4; up to 4 sets of distances can be measured. The results are displayed on the right-hand screen.


Fig. 1: Distance measurement

The 4 distances are D1, D2, D3 and D4.
D1/D2 is the ratio of D1 to D2. D3/D4 is the ratio of D3 to D4.

4.7 Volume measurement
There are two ways to measure the volume of a saved photo.

1. 3-set volume measurement

3-set spacer measurements are those in which the volume is calculated using a 3-axis method. The spacer must be measured three times before the volume can be calculated. Press M/0 to obtain the result.
If there are fewer than 3 sets of values, no result will be displayed when you press M/0. If you enter 4 sets and then press M/0, the result will comply with the first three values entered. Method: (see the image of a kidney below as an example):
1. Select the diagonal and the longitudinal section of the kidney and switch on Freeze mode.
2. Measure the length and short axes of the diagonal section using the distance measurement values.
3. Measure the diameter of the longitudinal section using the distance measurement function.
4. Press M/0 to complete the measurement with the 'Vm1' volume values on the right-hand side:


Fig. Volume measurement (3-set volume measurement)

2. 2-set volume measurement

Measurements using 2 sets of circumference and an area using the ellipse method. The result is compiled from the measurements. Method:
(see photo of a kidney below):
1. Take the diagonal and the longitudinal section of the kidney and switch on Freeze mode.
2. Measure the length and short axes of the diagonal section using the distance measurement values.
3. The system automatically completes the measurement, displaying the 'Vm1' volume values on the right-hand side.

4.8 Heartbeat measurement
1. Freeze the desired image in B/M mode
2. Measure the distance between the two largest waves over two periods using the distance measurement method.
3. A set of values appears in the bottom right-hand corner. The values are time T and heart rate HR.


Fig. 1: Heart rate measurement

4.8.1 Heart rate measurements using the spacer method
Use the left mouse button to activate the "+" cursor. Place the "+" at the starting point of the measurement, then use the right mouse button to fix the cursor in place. Place the second cursor at the end point of the measurement (if you right-click, you can switch between the start and end points). Press the left mouse button to measure the distance. Repeat the process to add further measurements.

4.8.2 Heart rate measurements using the circumference method
Use the left mouse button to activate the cursor. Move the cursor to the starting point of the area to be measured. Then click the middle mouse button to fix the cursor at the starting point. Move the cursor around the area to be measured. Use the middle mouse button to obtain the value. Repeat the process to add further measurements.

4.9 Estimating the gestation period in pigs
1. The pig must be standing for this.
2. Place the probe slightly to the left or right of centre on the ventral abdominal wall, running tightly alongside the teats and the skull down to the hind leg. If there is dirt at this point, clean it first with water if the abdominal-pelvic structure cannot be visualised correctly.




Fig. 1: Sketch showing the position during gestation


3. To measure the cardiac macroaxis, the screen should be large enough to display the maximum longitudinal axis of the heart. As gestation progresses, the macroaxis of the foetal heart increases steadily. The measurement method is shown in the following fig.:
4. Measure the selected parameter distance in accordance with the distance measurement method; the stage of pregnancy is automatically displayed after 'G.A.'.

Not pregnant

25 days after mating

90 days after insemination

4.10 Other

4.10.1 Printing
Connect the "Video Out" to the printer. (Cable sold separately). Switch on the device and follow the instructions

4.10.2 Transferring images to a PC
Connect one end of the USB cable to the device's USB port and the other end to the PC's USB port. Real-time images are transferred to the PC via the High-Speed USB 2.0 usb interface. For further information, please refer to the user manual on the CD.

4.10.3 Switching off the device
Do not plug in or unplug the mains plug; if the device needs to be restarted after being switched off, wait 2–3 minutes.
The device continues to consume power even in standby mode. If the device is not going to be used for a long period, please remove the battery from its compartment.

5. System maintenance
System maintenance is carried out by the user. Following purchase, the user is responsible for the operation and maintenance of the system. Check the condition of the device regularly and carry out preventative checks.

⚠ Warning: The following maintenance tasks must not be carried out during the examination.

5.1 Daily maintenance
Daily maintenance must be carried out by the user.

5.2 Cleaning and Disinfection

⚠ Warning: Switch off the appliance before cleaning it. Cleaning the appliance whilst it is switched on may cause damage to it.

⚠ Warning: The device has an IPX4 protection rating (no adverse effects from water splashes from any direction). Avoid splashing water or liquids onto the device, as this may cause damage.

⚠ important:
1. To prevent the possible transmission of diseases, it is recommended that you wear sterile gloves when performing cleaning and disinfecting tasks.
2. Avoid overheating the probe (above 60°C) during cleaning and disinfection, as this could cause it to become misshapen or damaged.

1. Cleaning the device
 (1) Wear sterile gloves to prevent possible contamination with pathogens.
 (2) Wipe the probe with water or soapy water to remove all dirt. Do not use brushes or other abrasive cleaning implements, as this may damage the probe.
 (3) After wiping, use a sterile cloth to wipe the water off the probe. Do not perform drying on the probe by heating or exposure to sunlight.
2. Disinfection of the device
 (1) Sterile gloves must be worn during routine disinfection to prevent any possible contamination.
 (2) Perform cleaning of the probe before routine disinfection, then use a 75% alcohol solution for repeated wiping of the probe.
 (3) Rinse any remaining chemicals off the probe with sterile water.
 (4) Use a sterile cloth to wipe the water from the surface of the probe. Do not perform drying on the probe by heating or exposure to sunlight.

⚠ important:
1. It is standard for the ultrasound focus to discolour and for the text on the device label to fade.
2. The number of routine disinfections should be kept to a minimum, as this can affect the safety and performance of the probe.

3. Inspection of the device after cleaning and disinfection
 (1) Check the device's casing for any irregularities (e.g. wear and tear, cracks or flaking).
 (2) Check the device's transducer for any irregularities such as wear, cracks, delamination, bulging and chipping.

5.3 Use for and replacement of the rechargeable battery

Use only batteries supplied by Kerbl. Do not remove the built-in rechargeable battery to avoid damage to the battery or the device. The battery should only be replaced by service personnel from Albert Kerbl GmbH.

5.4 Maintenance of the rechargeable battery

The battery is a consumable item; its service life depends on the number of charge and discharge cycles. If the operating time is significantly reduced compared with normal operation, the battery should be replaced immediately. Temperatures that are too high or too low impair charging and discharging performance and shorten the battery's service life and capacity.

⚠ Important: The charger must comply with CCC specifications or national requirements for chargers.

⚠ Important: Do not simply dispose of the faulty battery; in particular, do not throw it into a fire. Please dispose of the battery in accordance with local regulations.

important:

1. **Do not throw the battery into water or allow it to get wet, as this may cause leakage, explosions or fires.**
2. **Do not use or store the battery near heat sources such as fires or heaters, as this may cause leakage, explosions or fires.**
3. **Do not short-circuit the anode and cathode, as this may cause leakage, explosions or fires.**
4. **Do not heat the battery or throw it into a fire, as this may cause leakage, explosions or fires.**
5. **Do not connect the anode and cathode to metal or other conductors. Do not transport or store the battery together with neck chains, hairpins or other metallic objects, as this may cause leakage, explosions or fires.**
6. **Do not throw or shake the battery, as this may cause leakage, explosions or fires.**
7. **Do not insert nails or other sharp objects into the battery; do not step on the battery, as this may cause leakage, explosions or fires.**
8. **Do not clip the battery terminals together, as this may cause leakage, explosions or fires.**
9. **Do not dismantle the battery in any way, as this may cause leakage, explosions or fires.**
10. **Do not charge the battery near heat sources or in overheated environments, as this may cause leakage, explosions or fires.**
11. **Do not place the battery in a microwave or in a pressurised container, as this may cause leakage, explosions or fires.**
12. **Do not use the battery together with disposable batteries (such as dry batteries) or rechargeable batteries of different capacities, models or brands, as this may result in leakage, explosions or fires.**
13. **Do not use a defective battery that emits an unusual odour, generates excessive heat, is deformed, discoloured or exhibits any other abnormal characteristics, as this may result in leakage, explosions or fires.**
14. **Stop charging immediately and unplug the charging cable if the battery shows any abnormal signs, such as an unusual odour, excessive heat, deformation or discolouration. Failure to do so may result in leakage, explosions or fires.**
15. **Keep the battery away from sources of fire if leakage or unusual odours occur, as this may lead to leakage, explosions or fires.**
16. **If electrolyte splashes into your eye, do not rub it; rinse it out immediately and seek medical attention as soon as possible. Failure to do so may cause permanent damage to the eye.**
17. **Do not use the battery in hot environments, such as in direct sunlight or in a car at high temperatures, as this may cause fires, impair output and shorten its service life.**
18. **If you do not use the battery under the conditions specified in this manual, this may impair output, shorten its service life or even lead to extreme temperatures, explosions or fires.**

6. General instructions

6.1 Warranty and Repairs

Albert Kerbl GmbH provides a 1-year warranty on this item in accordance with our general warranty terms and conditions. These general warranty terms and conditions can be found at www.kerbl.com/guarantee. The warranty and statutory guarantee will lapse in the event of improper use. Damaged parts must only be repaired by a qualified technician or the manufacturer. Only use original spare parts for any repairs.

6.2 Electronic Scrap

Do not dispose of the battery in household waste. Dispose of the battery properly, separately from the electrical appliance and in accordance with local regulations. Municipal collection points and take-back points at retailers are available for this purpose. The separate collection and recycling of old batteries and accumulators is intended to prevent adverse effects on the environment and human health.

6.3 CE Declaration of Conformity

CE

Albert KERBL GmbHhereby declares that the product/device described in this manual complies with the essential requirements and all other relevant national and European regulations. The CE mark indicates compliance with European Union directives. The Declaration of Conformity can be viewed at the following web address: www.kerbl.com.

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Kerbl UK Ltd 8 Lands End Way