Ultraschallscanner EXPERT inkl. Rektalsonde #221302 Bedienungsanleitung

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2. Installation 24
2.1 Ambient conditions 24
2.2 General instructions for operating the device 24
2.3 Installation and removal of the probe 25

3. Operation 25
3.1 Display 25
3.2 Functions of the control buttons 26

4. Operating instructions 26
4.1 Switching on 26
4.2 Diagnostic procedure 26
4.3 Freezing the image 26
4.4 Adjusting image parameters 26
4.4.1 Gain setting 26
4.4.2 Brightness 26
4.4.3 Contrast 26
4.4.4 Image smoothing 26
4.4.5 Frequency 27
4.4.6 Energy-saving mode 27
4.5 Function descriptions 27
4.5.1 Cine-Loop image playback 27
4.5.2 Image storage 27
4.5.3 Image playback 27
4.5.4 Deleting images 27
4.5.5 Backfat/Length Measurement 27
4.5.6 Reading from USB storage media 27
4.5.7 Operating modes 27
4.5.8 Gamma correction 27
4.6 Performing examinations 28
4.6.1 Examination of horses (rectal probe) 28
4.6.2 Examination of cattle (rectal probe) 28
4.6.3 Examination in pigs (convex probe) 28
4.6.4 Examination in sheep (convex probe) 29

5. System maintenance 29
5.1 Daily maintenance 29
5.2 Cleaning and Disinfection 29
1. Cleaning the device 29
2. Disinfection of the device 29
3. Inspection of the device after cleaning and disinfection 30
5.3 Use only batteries supplied by Kerbl. 30
5.4 Maintenance of the rechargeable battery 30

6. Troubleshooting and Solutions 31

7. General instructions 31
7.1 Warranty and repairs 31
7.2 Electronic Scrap 31
7.3 Declaration of Conformity 31

(EN) Operating instructions
EXPERT Ultrasound Scanner including Rectal Probe #221302

1. Product overview

1.1 Intended use

The EXPERT ultrasound scanner is a portable ultrasound device for veterinary applications. It consists of a main unit and can be used with a wired 6.5 MHz electronic linear rectal probe for ultrasound examinations of cattle or horses, or with a 3.5 MHz convex probe for ultrasound examinations of pigs, sheep, goats, dogs, cats or other pets.

1.2 Device overview

Front view

Side view with interfaces

1.3 Technical Data

	221302-2 6.5 MHz rectal linear probe	221302-1 3.5 MHz convex probe
Scan depth (mm)	≥80	≥120
Transducer	96 elements	96 elements
Resolution (mm) – lateral	≤1 (depth ≤60 mm)	≤3 (depth ≤80 mm), ≤5 (80 mm < depth ≤130 mm)
Resolution (mm) – axial	≤1 (depth ≤80 mm)	≤1 (depth ≤80 mm)
Blind zone (mm)	≤3	≤6
Geometric positional accuracy – horizontal	≤5 %	≤7.5 %
Geometric positional accuracy – vertical	≤5 %	≤5 %
Screen size	8-inch LCD	
Resolution	1024x720 pixels	
Display modes	B, B+B, 4B, B+M, M	
Grey scales	256	
Permanent image memory	8 GB	
Gain range	0–200 dB (adjustable range 64–200 dB)	
Cine-Loop memory	≥500 individual frames	
Adjustable scan depth	80–240 mm	
Gain control	Overall gain, near-field, mid-field and far-field gain	
Focus control	Number of focus points, focus position	
Image processing	Pseudo-colours, greyscale calibration, image smoothing, histogram	
Adjustable frequencies	4 / 5 / 6.5 / 7.5 / 8 MHz	2.5 / 3 / 3.5 / 4 / 4.5 / 5 / 5.5 MHz
Frame correlation	Supported	
Measurement functions	Length, circumference, area, volume	
Connections	USB-A	
Water resistance	IP X6	
Battery life	> 4 hours	
Device dimensions	225 × 213 × 60 mm	
Device weight	1229 g	
power consumption	9 W	
battery	18650 / 11.1 V / 6800 mAh	

1.4 Basic principle

B-mode ultrasound provides visibility into internal structures by sending ultrasound waves into the tissue and analysing the returning echoes. This works as follows: Different types of tissue have different densities and therefore different levels of resistance to ultrasound, i.e. different acoustic impedances (a combination of sound density and sound velocity). When the piezoelectric chip (transducer) receives a specific electrical impulse, it generates ultrasound at a specific frequency. When these ultrasound waves penetrate the tissue, the various organs produce different echoes. The differently reflected sound waves are received again by the transducer, which converts them into electrical impulses, reinforces them, organises them and digitises them. A standard video signal is generated and cross-sectional images of organs are displayed on the screen.

1.5 Standard configuration

221302-8	EXPERT main unit
221302-2	6.5 MHz rectal linear probe
221302-7	6800 mAh Li-ion battery, 11.1V
	Neck strap with 4 carabiners
221302-5	Carry case
221301-7	Charger (EU) for PRO & EXPERT
221301-8	UK plug with cable for PRO & EXPERT

1.6 Optional accessories

221302-1	3.5 MHz convex probe
221302-6	Anti-glare screen with collar
221302-4	Extension handle for rectal probe

⚠ Warning: Please use only the original parts specified above. Non-original parts may cause unforeseen electromagnetic compatibility (EMC) issues. Any potential EMC coupling effects must be approved separately

2. Installation

2.1 Ambient conditions

1. Ambient temperature range: +5 °C to +40 °C
2. Ambient humidity: 30% to 75%
3. Atmospheric pressure range: 70 kPa to 106 kPa
4. Mains voltage: 100–240 V AC, 50 Hz ±1 Hz / 60 Hz ±1 Hz

Avoid severe vibrations during operation and keep the device away from equipment generating strong electric fields, strong magnetic fields or high voltage. Do not expose the display to direct sunlight. Ensure the device is well ventilated and protect it from moisture and dust.
After unpacking, check the 'Packing List' for any transport damage, then proceed with the installation.

⚠ Warning:
1. If the device shows signs of damage caused during transport after unpacking, do not put it into operation.
2. Handle the probes with care. Damage caused by dropping or mechanical pressure is not covered by the manufacturer's warranty.

2.2 General instructions for operating the device

During operation

1. Do not switch the device back on within 2–3 minutes of switching it off.
2. If you notice anything unusual during the examination, stop the measurement immediately and switch off the device.

After use


1. Switch off the device.
2. Pull out the power plug from the socket – do not pull on the cable.
3. Remove the ultrasound gel from the probe using soft, sterile medical cotton wool pads.

2.3 Installation and removal of the probe

There is only one probe socket to which all supported optional probes can be connected.

Installation of the probe
Align the recess on the probe connector with the socket on the main unit, plug the connector in and tighten the metal housing.
Dismantling is carried out in the reverse order to installation.


Installation of the probe


⚠ Warning: Under no circumstancesshould you disconnect or connect the probe plug whilst the device is switched on if the probe or the main unit is damaged. Once the probe is connected to the main unit, do not simply disconnect or reconnect it, even if there are connection problems.

3. Operation

3.1 Display

3.2 Functions of the control buttons

4. Operating instructions

4.1 Switching on

Press the ' ' to start the device. Then press any button to access the system.

Note: The display quality of the LCD screen depends heavily on the viewing angle; adjust this if necessary.

4.2 Diagnostic procedure

Apply a sufficient amount of ultrasound gel to the area to be examined on the animal. Press the probe's transducer flat against the area. A cross-sectional image of the tissue will appear on the screen. Move the probe carefully to set the optimum scanning depth. At the same time, adjust the overall gain, near-field/far-field gain or the focus to obtain the best possible image.

Note

1. Do not press the probe too hard against the animal – this may damage the probe or cause the animal pain.

2. Always use the probe with the correct frequency suitable for the area being examined.

4.3 Freezing the image

Press ' ' to switch between the live image and a frozen image.

Note: The message 'FREEZE' appears below in the right-hand corner of the screen when the image is frozen.

4.4 Adjusting image parameters

Adjustable parameters: probe frequency, gain, brightness, image smoothing, frame, scan angle and focus.

In real-time mode, press the " " button – a menu item will be highlighted. Use the " " and " " buttons to adjust the values; the settings are displayed on the right-hand side of the screen.

4.4.1 Gain setting

1. In real-time mode, press the menu button and use the ' '/' ' buttons to select the overall gain, near-field, medium-field or far-field gain (displayed in the top right-hand corner of the screen).
 Use the ' '/' ' buttons to change the values.
 - Overall gain: 0–100
 - Near-field/far-field gain: 0–30
2. Press " " twice to exit.

4.4.2 Brightness

1. In real-time mode, press " ", then use " "/" " to select "Brightness". Use the " "/" " buttons to adjust (range 0–63).
2. Press " " twice to exit.

4.4.3 Contrast

1. In real-time mode, press " ", then use the " "/" " buttons to select "Contrast". Use the " "/" " buttons to adjust (range 0–63).
2. Press " " twice to exit.

4.4.4 Image smoothing

1. Press " " again, then use the " "/" " buttons to select "Image smoothing". Use the " "/" " buttons to adjust the setting (range 0–4).
2. Press " " twice to exit.

4.4.5 Frequency

1. In real-time mode, press " ", then use the " "/" " buttons to select "Frequency". Adjust using the " "/" " buttons.
 - 3.5 MHz convex probe: 2.5 / 3 / 3.5 / 4 / 4.5 / 5 / 5.5 MHz
 - 6.5 MHz rectal probe: 4 / 5 / 6.5 / 7.5 / 8 MHz
2. Press ' ' twice to exit.

4.4.6 Energy-saving mode

1. In real-time mode, press " ", then use " "/" " to select "Energy-saving mode". Use " "/" " to set the time. Four settings available: off, 5 min., 10 min., 15 min.
2. Press " " twice to exit.

4.5 Function descriptions

4.5.1 Cine-Loop image playback

Images are temporarily cached (lost when the battery is empty) and then played back in an endless loop.
Procedure:

1. In real-time mode, recording for the Cine-Loop starts automatically; the recording lasts 30 seconds.
2. Press the following buttons in sequence: " "→" " (to activate "Cine Loop") → " " to start playback. The status "120/500" is displayed below in the left-hand corner (current frame / total number of frames stored).
3. Press the " " or " " buttons to activate single-frame playback.
4. Press the " " button to stop the film loop.

Note

If you change the scan method, probe or display mode, it will take another 30 seconds for the Cine-Loop to start recording again. Changes to the scan angle also affect the recording duration; the status is displayed in the bottom left-hand corner.

4.5.2 Image storage

The device has 8 GB of non-volatile memory (images are retained even in the event of a power failure).
Procedure for saving:

1. Carry out the examination and freeze a suitable image.
2. " "→ " " (select "Save & Display") → " "→Select "Save" → " ". "SAVING..." appears in the top left-hand corner – wait until "OK" is displayed in the bottom left-hand corner before making any further entries.
3. Press the " " to return to live view mode.

4.5.3 Image playback

1. " " → " " (select "Save & View") → " " → Select "View" → " ".
2. Use " "/" " to switch between saved images.
3. Press " " to return to live view mode.

4.5.4 Deleting images

1. " " → " " (select "Save & View") → " " → Select "Delete" or "Delete All" → " ", select the images to be deleted.
2. Press the ' ' button to return to real-time mode.

4.5.5 Backfat/Length Measurement

Backfat measurement

1. Press " " →, select "Measurement" →, " ".
2. Use the arrow keys to position the cursor; confirm the start and end points with " ".
3. The measured backfat value is displayed on the right-hand side of the screen.

Length, circumference, area and volume measurement

- When measuring distance, press the " " button once to switch to circumference and area measurement.
- Press the button again: volume measurement

Special animal measurements for determining pregnancy:

A.CRL (crown-rump length): Accurate pregnancy assessment in the early stages of pregnancy
B.BPD (biparietal diameter): Gestation determination in the medium stage

4.5.6 Reading from USB storage media

In real-time mode, connect the device to the computer using a dual-ended adaptor cable that connects both sides. As soon as the computer recognises the USB drive, the user can access the images stored on it: Computer → SanDisk → USB drive

4.5.7 Operating modes

Press ' ' twice in real-time mode – imaging modes are available for selection:

1. B-mode: standard ultrasound imaging
2. B\B: Dual B-scan split screen
3. B\M: B-mode + M-motion curve
4. M-mode: Motion curve only

4.5.8 Gamma correction


Press ' ' twice in real-time mode, select 'Gamma correction' using ' '/' ' and press ' '. Four profiles available:

1. Standard gamma: Balanced light-dark display, with universal use in pregnancy and organ examinations.
2. High Gamma: Dark areas emphasised, bright areas compressed – ideal for deep lesions and homogeneous tissues.
3. Extremely high gamma: Significant enhancement of dark areas, significant attenuation of bright areas – for hypoechoic structures (fluid collections, cysts).
4. Low gamma: Bright areas emphasised, dark areas compressed – a brighter overall image for superficial, high-echo structures (shadows are lost).

4.6 Performing examinations

4.6.1 Examination of horses (rectal probe)


1. Remove faeces from the horse's rectum.
2. Perform manual palpation to make a coarse assessment of pregnancy status in order to locate the tissue areas to be examined.
3. Grasp the probe and carefully insert it into the rectum; place your hand between the rectal wall and the probe head to immediately detect any unusual sensations.
4. The internal organs are displayed on the screen: the bladder is located at the vertical junction, followed by the cervix and the body of the uterus. The round uterine horns are visible horizontally. Move the probe within the rectum until the transition between the body of the uterus and the uterine horn is visible, then rotate the probe towards the uterine horn.

Schematic representation of the probe position in the uterus and ovary
1 Rectum | 2 Uterine horn | 3 Body of the uterus | 4 Ovary | 5 Vagina | 6 Urinary bladder

4.6.2 Examination of cattle (rectal probe)


1. Remove faeces from the cow's rectum.
2. Perform manual palpation to make a coarse assessment of pregnancy status.
3. Carefully insert the probe into the rectum, placing your hand between the rectal wall and the probe.
4. Organs are displayed on the screen: bladder, cervix, body of the uterus and round uterine horns. Move the probe until the transition between the body of the uterus and the uterine horn is visible; rotate the probe towards the uterine horn.


Schematic representation of the probe position in the uterus and ovary
1 Rectum | 2 Uterine horn | 3 Body of the uterus | 4 Ovary | 5 Vagina | 6 Urinary bladder

4.6.3 Examination in pigs (convex probe)


1. Keep the pig in a standing position; perform cleaning of any heavily soiled skin at the examination site beforehand to ensure better image quality.
2. Optimal examination site: 10–15 cm above the second pair of teats on the abdomen. Place the probe on the uterine area, tilted towards the front at a 45° angle.
3. Black, round areas on the image indicate either incorrect probe positioning or a full bladder – allow the pig to urinate and repeat the examination.
4. Not pregnant: The uterus appears white.
 After mating: Early pregnancy shows black amniotic sacs; late pregnancy shows white bony structures.


Schematic examination position in pigs

Not pregnant

25 days after mating

90 days after insemination

4.6.4 Examination in sheep (convex probe)

Two common B-ultrasound methods for pregnancy diagnosis: abdominal examination with a convex/linear probe or rectal examination with a rectal probe. Both methods are reliable:

- Up to 35 days' gestation: Rectal examination is more accurate than abdominal examination
- 35–70 days' gestation: Both methods are equally effective
- Over 70 days' gestation: Abdominal examination is more practical due to the enlarged uterus

Abdominal examination

1. Examine the sheep whilst standing, lying down or crouching. Place the probe on the hairless central abdominal area and position it correctly.
2. Remove any mud or dirt from the area with water to ensure a clear view of the pelvic organs.

5. System maintenance

System maintenance is carried out by the user. Following purchase, the user is responsible for the operation and maintenance of the system. Check the condition of the device regularly and carry out preventative checks.

⚠ Warning:
The following maintenance tasks must not be carried out during the inspection

5.1 Daily maintenance

Daily maintenance must be carried out by the user

5.2 Cleaning and Disinfection

⚠ Warning:
Switch off the appliance before cleaning it. Cleaning the appliance whilst it is switched on may cause damage to the appliance.

⚠ Warning: The device has an IPX6 protection rating (no adverse effects from water spray from any direction). Avoid spraying water or liquids onto the device, as this may cause damage.

⚠ important:

1. **To prevent the possible transmission of diseases, it is recommended that you wear sterile gloves when performing cleaning and disinfecting.**
2. **Avoid overheating the probe (above 60°C) during cleaning and disinfection, as this could cause it to become misshapen or damaged.**

1. Cleaning the device

1. Wear sterile gloves to prevent possible contamination with pathogens.
2. Wipe the probe with water or soapy water to remove all dirt. Do not use brushes or other abrasive cleaning implements, as this may damage the probe.
3. After wiping, use a sterile cloth to wipe the water off the probe. Do not perform drying on the probe by heating or exposure to sunlight.

2. Disinfection of the device

1. Sterile gloves must be worn during routine disinfection to prevent any potential contamination.
2. Perform cleaning of the probe before routine disinfection, then use a 75% alcohol solution for repeated wiping of the probe.
3. Rinse any remaining chemicals off the probe with sterile water.
4. Use a sterile cloth to wipe the water from the surface of the probe. Do not perform drying on the probe using heat or sunlight.

⚠ Important:

1. **It is standard for the ultrasound focus to discolour and for the text on the device label to fade.**
2. **The number of routine disinfections should be kept to a minimum, as this can affect the safety and output of the probe.**

3. Inspection of the device after cleaning and disinfection

1. Check the device's casing for any irregularities (e.g. wear and tear, cracks or flaking).
2. Check the device's transducer for any irregularities such as wear, cracks, delamination, bulging and chipping.

5.3 Use only batteries supplied by Kerbl.

5.4 Maintenance of the rechargeable battery

5.4.1 The battery is a consumable item; its service life depends on the number of charge and discharge cycles. If the operating time is significantly reduced compared with standard operation, the battery should be replaced immediately.

5.4.2 Temperatures that are too high or too low impair charging and discharging performance and reduce the battery's service life and capacity

⚠ Important: The charger must comply with CCC specifications or national requirements for chargers.

⚠ Important: Do not simply dispose of the faulty battery; in particular, do not throw it into a fire. Please dispose of the battery in accordance with local regulations.

⚠ important:

1. **Do not throw the battery into water or allow it to get wet, as this may cause leakage, explosions or fires.**
2. **Do not use or store the battery near heat sources such as fires or heaters, as this may cause leakage, explosions or fires.**
3. **Do not short the anode and cathode, as this may cause leakage, explosions or fires.**
4. **Do not heat the battery or throw it into a fire, as this may cause leakage, explosions or fires.**
5. **Do not connect the anode and cathode to metal or other conductors. Do not transport or store the battery together with neck chains, hairpins or other metallic objects, as this may result in leakage, explosions or fires.**
6. **Do not throw or shake the battery, as this may cause leakage, explosions or fires.**
7. **Do not insert nails or other sharp objects into the battery; do not step on the battery, as this may cause leakage, explosions or fires.**
8. **Do not clip the battery terminals together, as this may cause leakage, explosions or fires.**
9. **Do not dismantle the battery in any way, as this may result in leakage, explosions or fires.**
10. **Do not charge the battery near heat sources or in overheated environments, as this may cause leakage, explosions or fires.**
11. **Do not place the battery in a microwave or in a pressurised container, as this may cause leakage, explosions or fires.**
12. **Do not use the battery together with disposable batteries (such as dry batteries) or rechargeable batteries of different capacities, models or brands, as this may result in leakage, explosions or fires.**
13. **Do not use a defective battery that emits an unusual odour, generates excessive heat, is deformed, discoloured or exhibits any other abnormal characteristics, as this may result in leakage, explosions or fires.**
14. **Stop charging immediately and unplug the charging cable if the battery shows any unusual signs, such as an unusual odour, excessive heat, deformation or discolouration. Failure to do so may result in leakage, explosions or fires.**
15. **Keep the battery away from sources of fire if leakage or unusual odours occur, as this may lead to leakage, explosions or fires.**
16. **If electrolyte splashes into your eye, do not rub it; rinse it out immediately and seek medical attention as soon as possible. Failure to do so may cause permanent damage to the eye.**
17. **Do not use the battery in hot environments, such as in direct sunlight or in a car at high temperatures, as this may cause fires, impair output and shorten the battery's service life.**
18. **Failure to use the battery under the conditions specified in this manual may impair output, shorten its service life, or even lead to extreme temperatures, explosions or fires.**

6. Troubleshooting and Solutions

- Simple faults can be rectified simply by the user using the table below. For complex faults, contact the manufacturer.
- The device must only be opened and repaired by trained specialist personnel.

No.	Fault description	Suggested solutions
1	Power supply indicator is not lit	Check the mains supply Check the mains cable and socket
2	Device power indicator is off	Check the device's power cable and socket
3	Screen has streaks or a snowy, distorted image	Check the power supply Check for electromagnetic interference from the surroundings Check the probe socket
4	Blurry image	Adjust the overall gain Adjust brightness and contrast
5	Delayed image	Press the '⊜' menu button during operation Press the Up/Down button to navigate to 'Frame C' Use the Left/Right buttons to optimally adjust 'Frame C': → The lower the value, the shorter the image processing time and the lower the latency

7. General instructions

7.1 Warranty and repairs

Albert Kerbl GmbH provides a 1-year warranty on this item in accordance with our general warranty terms and conditions. These general warranty terms and conditions can be found at www.kerbl.com/guarantee. The warranty and statutory guarantee will lapse in the event of improper use. Have damaged parts repaired only by a qualified technician or the manufacturer. Use only original spare parts for any repairs.

7.2 Electronic Scrap

Do not dispose of the battery in household waste. Dispose of the battery properly, separately from the electrical appliance and in accordance with local regulations. Municipal collection points and take-back points at retailers are available for this purpose. The separate collection and recycling of old batteries and accumulators is intended to prevent adverse effects on the environment and human health.

7.3 Declaration of Conformity

Albert KERBL GmbH hereby declares that the product/device described in this manual complies with the essential requirements and all other relevant national and European regulations. The CE mark indicates compliance with European Union directives. The Declaration of Conformity can be viewed at the following web address: www.kerbl.com.

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Kerbl UK Ltd 8 Lands End Way