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(EN) Oestrus detection

Product overview

Successful reproduction management relies heavily on the reliable identification of the optimal insemination time. The Rheintechnik brand's electronic oestrus detection helps farmers to determine oestrus more reliably and to narrow down the suitable time for insemination.

The device is particularly suitable for assisting with heat detection in cattle and cows with oestrus symptoms that are difficult to recognise (silent oestrus). Regular measurements enable changes in the course of oestrus to be detected and documented at an early stage.

The oestrus detection measures the electrical resistance of the vaginal mucus. This changes during the oestrus cycle depending on the animal's hormonal balance.

To ensure a reliable evaluation, measurements should always be carried out under comparable conditions and at regular intervals. The measurement results must always be evaluated in conjunction with observations of the animal's behaviour and other signs of oestrus.

Commissioning

Do not use the device if there is any visible damage.

1. When the device is switched on, two flashing lines appear.
2. After pressing the switch three times, the numbers 1, 2, 3 appear on the screen, indicating that three consecutive measurements have been recorded.
3. After pressing the switch again, the device displays a result that complies with the result of a '1 0' measurement cycle, indicating that the battery level is sufficient, the device is functioning correctly and the probe is not immersed in a solution. (i.e. the measurements were taken in air)

The smallest unit displayed is 10; the total measuring range is between 0 and 1990 units. If the high value is exceeded, the numbers '1 0' are displayed, similar to a measurement taken in air. The device's measuring range is many times larger than the maximum resistance value of vaginal mucus. If the measurement is carried out correctly, it is not possible to exceed the measuring range.

If necessary, you can check that the device is working properly using saline solution as follows:

1. Immerse the probe vertically into a glass of clean water at room temperature. Allow the probe to rest in the centre of the glass base, carry out a measurement cycle and read the result. It will probably display '1 0', which means that the measuring range has been exceeded, as the resistance of clean water is very high.
2. Next, add a small pinch of salt, mix it well into the water, immerse the probe as described in point v, and carry out the measurement cycle again. After adding the salt, the result will change depending on the salt concentration, as salt reduces resistance (lower resistance = lower readings).

Measurements on animals


1. Before you start the measurement, ensure that the electrodes are clean and free from grease.
2. You should not use any lubricants whilst using the oestrus detection system, as they interfere with the electrodes' contact with vaginal mucus and distort the final result.
3. Is the battery charged?
4. Prepare a disinfectant solution in accordance with the manufacturer's instructions to sterilise the probe before and after use (see the 'Disinfection' section for further details)
5. If the vulvar area is soiled, perform cleaning to ensure accurate measurement results.
6. If necessary, ask another person to assist you with the measurement.

Spread the animal's vulva and carefully begin to insert the disinfected probe into the vagina.

It is best to insert the probe from below to above at a 45-degree angle to prevent the probe from entering the urinary tract.

Insert the probe until you feel slight resistance. At this point, the probe is positioned at the cervix. Then make two or three half-turns to ensure the best possible contact between the electrodes and the oestrus mucus.

The measurement should always be carried out in the same way.

NOTE! The insertion depth of the probe depends on the length of the vagina and may vary between animals; you should therefore always insert the probe until you feel slight resistance.

1. Press the on/off switch. As with the test measurement, two flashing lines will be displayed.

2. Carry out a measurement cycle and remember to make 2–3 half-turns between each measurement, each time after pressing the power button.
3. An example of a measurement cycle:

4. Remove the probe from the animal's vagina and record the result.
5. Switch off the device by holding down the on/off switch for more than 3 seconds, or wait until it switches off automatically after 30 seconds.
6. Wash and disinfect the device as described in the description of disinfection.

Heat detection


Fig. 1 shows the change in the resistance of vaginal mucus over the course of the ovarian cycle. During the non-oestrus period, a high resistance value is measured (depending on the type of animal, e.g. approx. 300 units or more in cows).
As ovulation approaches, the resistance decreases and reaches its minimum value (e.g. around 200 units in cows). From ovulation onwards, it rises again to a high level and remains at this level until the next ovulation.




Fig. 2 shows the progression of resistance following successful insemination. The resistance value remains at a consistently high level throughout pregnancy until the due date.

The following diagrams illustrate how the electrical resistance of vaginal mucus changes during the oestrus cycle in cows and mares. For heat detection, the lowest reading and the subsequent marked increase in resistance are particularly crucial. This increase indicates impending ovulation and helps to determine the optimal time for insemination or mating. Monitor changes in the readings over several days and always carry out the measurements under comparable conditions. As the oestrus phase is short-lived, it is advisable to take measurements several times a day during this period.
Note: The resistance values shown in the diagrams are for guidance only. The actual readings may vary depending on the animal, its age, state of health and housing conditions. It is therefore the trend in the readings that is crucial, rather than the absolute numerical value.


Disinfection

Thorough cleaning and disinfection are essential for the reliable and hygienic use of the oestrus detection system.
1. Disinfect the oestrus detection before and after each measurement.
2. Dirty or greasy electrodes can distort the measurement results and lead to elevated readings.
3. Urine or salty liquids can also affect the measurement results and lead to lower readings.

Clean the probe thoroughly with water after use and then immerse it in a suitable disinfectant solution.
Note: Do not perform cleaning with very hot or boiling water.

Prepare the disinfectant solution in a suitable container. Use two separate containers: one for cleaning and one for disinfecting the device. The concentration of the disinfectant must comply with the manufacturer's instructions. Solutions with a concentration of 0.5–1 per cent are often recommended. Immerse the probe in the solution for several seconds to ensure the disinfectant is reliable.

Note: Disinfectants can irritate the animals' mucous membranes. It is therefore essential to adhere to the recommended concentration. With frequent use, the solution loses its effectiveness over time and should be freshly prepared regularly. Poor hygiene or the incorrect use of disinfectants can cause irritation, inflammation or damage to the animals' genital area.

Instructions for use

1. Use the device exclusively on your own farm for hygiene reasons. Passing it on to other farms may increase the risk of disease transmission if it is not cleaned or disinfected properly.
2. Always keep the oestrus detection system clean and disinfect it thoroughly before and after each application.
3. Store the device in a dry place at room temperature.
4. Do not use any lubricants for the measurement. These can interfere with the contact between the electrodes and the vaginal mucus and distort the measurement results.

5. Persistently low resistance readings over several days may indicate inflammation in the animals' genital area.
6. A damaged or missing protective cap may allow water to enter the device and cause permanent damage.
7. Do not open or remove the display yourself. This may damage the device's seal.

Changing the battery

1. Loosen the screws that fasten the flap to the handle.
2. Remove the battery from the battery compartment
3. Insert the new battery
4. When closing the battery compartment, ensure that the seal is correctly positioned.

Technical data

Weight	428 g incl. battery
Probe length	47 cm
Power supply	1 x 9 V 6F22 battery
Current consumption	Approx. 18 mA
Display	2.0" digital LCD display
Measuring range	0 – 1990s units
Operating and storage temperature	0° – +50°C
Maximum storage humidity	85 %
Housing material	ABS

Disposal

The operator is responsible for the proper disposal of the device once it is no longer in working order. Please observe the relevant regulations in your country. The device must not be disposed of with household waste. Under the EU Directive on the disposal of electrical and electronic equipment, the device accepted free of charge at local collection points or recycling centres, or it can be returned to specialist retailers offering a take-back service. Proper disposal helps protect the environment and prevents potential harmful effects on persons and the environment.

Declaration of conformity

Albert KERBL GmbH hereby declares that the product/device described in this manual complies with the essential requirements and all other relevant national and European regulations. The CE mark indicates compliance with European Union directives. The declaration of conformity can be viewed at the following web address: www.kerbl.com.

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Albert Kerbl GmbH
Kerbl UK Ltd Lands End Way